X SERIES® ADVANCED

Because Children Are Not Small Adults

Too often, medics have no choice but to adapt their equipment to treat sick or injured children. With the X Series® Advanced monitor/defibrillator, medics can optimize their treatment of young patients with specific modes for pediatrics and neonates.



Distinct Pediatric and Neonatal Modes

Children require specialized monitoring and defibrillation capabilities. The X Series Advanced automatically lowers the energy level of the shock and activates a pediatric-specific arrhythmia analysis algorithm simply by connecting to ZOLL® pediatric electrodes.

Pediatric CPR Feedback

To maximize CPR performance, the X Series Advanced CPR Dashboard™ displays the actual depth and rate of compression delivery in real time.

ZOLL Medical Corporation | 269 Mill Road | Chelmsford, MA 01824 | For contact details, visit www.zoll.com

© 2021 ZOLL Medical Corporation. All rights reserved. X Series, CPR Dashboard, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

